

# X-Ray Lab & Imaging Currents

Quarterly Publication of the American Association of Dental Maxillofacial Radiographic Technicians

Winter 2004

## Cone Beam Computed Tomography in Dentomaxillofacial Imaging

By: *Predrag Sukovic*

Various imaging modalities have been used in the dentomaxillofacial fields over the past few decades, none of them with entirely satisfactory results. This is particularly true for more demanding imaging tasks, such as implant planning, temporomandibular joint imaging, detection of facial fractures, lesions and diseases of soft tissue in the head and neck, and reconstructive facial surgery.

In particular, the use of dental implants is becoming the treatment of choice for the replacement of missing teeth. The successful outcome of a dental implant—the osseointegration of the implant—is heavily dependent on precise presurgical planning. Since the functional load on implants can be high, it is important that the implant be placed in a position where it can contact cortical bone and at an angle where the forces are as perpendicular as possible. Selection of the appropriate size and inclination of the implant in both a bucco-lingual and mesio-distal direction requires precise knowledge of the anatomy of the proposed site, including its dimension in all planes, the presence of knife-edge ridges and undercuts, as well as the location of anatomic structures, such as the nasal fossae, the maxillary sinus, and the mandibular canal. An evaluation of the thickness of the cortical bone and the density of the medullary bone is also critical to the success of the implant.

Commonly used dentomaxillofacial imaging modalities, such as periapical radiography, panoramic radiography, and conventional tomography produce only two dimensional and/or distorted images. As a result, a number of practitioners have resorted to outsourcing CT scans for implant planning and other demanding imaging tasks.

*Cone Beam Continued on page 10*



**Predrag Sukovic**  
President Xoran Technologies  
Ann Arbor Michigan, USA

## Inside This Issue

<b>Cone Beam Imaging</b>	<i>Predrag Sukovic</i> .....1
<b>President's Message</b>	<i>Craig Dial</i> .....2
<b>Editor's View</b>	<i>Camille Mayorga</i> .....5
<b>News and Trends</b>	.....6
<b>Spotlight Lab</b>	<i>Betty Walls</i> .....8
<b>CBCT Devices</b>	<i>CDA Journal</i> ..... 9
<b>AADMRT Committee Reports</b>	.....13
<b>Lab Products</b>	.....18

# President's Message



*Craig Dial*

It is only the second week of February, but here in California the weather is already beginning to change toward spring. I sometimes walk around the neighborhood of the lab during my lunch hour, and just this week I have begun to notice the changes this warmer weather brings. As I walk, I see a multitude of trees beginning to bud out, with their swelled-up tips ready to burst out a display of flowers and leaves. As I look down I notice plants and shrubs that are coming out of their deep sleep, ready to welcome the sun. Although it is technically winter, I know that many people around me can feel, smell and taste the spring in the air, and this has had an invigorating effect on me, and the people around me.

When we experience this phenomenon, we call it spring fever; with delight it comes over us, awakening our senses with rejuvenation. We feel good to be alive, and it excites our mind and body to move onward with our current projects, future plans, and projected goals. Just like the plants and animals that welcome this time of year, we too welcome this change, because change motivates us toward our ambitions and dreams.

Your part of the county (or world) may not be experiencing spring fever in the air, trees, or in your colleagues yet, but you can carry that rejuvenating spring in your step, in your voice, and in your heart. Think about coming out of this winter hibernation today, and commit yourself into an innovative view of life. You will be glad you did.



In front of DDI Sacramento (above), across the street (below)



## 2003 OFFICERS

American Association of Dental Maxillofacial Radiographic Technicians

### Executive Committee

President	Craig Dial
Vice President	Duane Perry
Secretary	Merry Hampton
Treasurer	Claudette Buehler

### Board Members

Membership	Kathleen Cox
Recruitment	Kevin Fox
Nominating	Jeannie Herriott
Advertising	Bart Webb

### Standing Committee Chairs

By Laws	Tracey Saucier
Continuing Education	Matt Kroona
Staff Photographer	Phil Abel
Research	Craig Dial & Don Croall
Web Master	Matt Kroona
Directory	Devery Wallace
Editor	Camille Mayorga



**Sharper.**

**We make it a point to give you more detail.** That's why Kodak developed two of its sharpest-resolution tools for dental practices yet—KODAK EKTAVISION G and L Extraoral Films. Whether it's our high-contrast G film, or our wide-latitude L film, you'll get low-crossover technology that captures

crisp radiographic detail, ideal for panoramic and cephalometric radiographs, TMJ cases, implant procedures, and more. Add in a 400-speed film system that minimizes radiation exposure, and you'll get an indispensable tool designed especially for your dental practice.

**But what else would you expect from the most trusted name in imaging?**

DENTAL



Following in the CommCAT tradition...

# PANOREX CMT and PANOREX CMT PLUS

Complex Motion Tomography for Implants,  
TMJ, Sinuses & more

Panoramic Imaging

Cephalometric Imaging



**Serving the AADMRT for over 10 years**

*Manufactured in the USA*



*Pending 510K approval*

# Cone Beam CT Technology



Full 3D Volumetric Scan  
with patient scan time of less than 40 seconds

Compact Design with the patient in a sitting position

Priced Affordably



800-205-3570, 215-997-5666  
[www.imagingsciences.com](http://www.imagingsciences.com)

## We are proud to introduce...



### PAN & CEPH DENTAL FILM

**Experience The Quality...Enjoy The Value**

# **Panoramic Dental U.S. Corp**

**Toll Free: 1-800-624-2345 · Toll Free Fax: 1-888-437-0200**  
**email: [lordjohn@shawcable.com](mailto:lordjohn@shawcable.com)**

# Editors View

With a new year ahead of us, we have the opportunity to refresh our ideas and weed out any habits that are not serving us. Now is a good time to re-evaluate the way we think about our lives, our jobs, and the way we perceive the world. It is a great time to let go of any negative thought patterns we have incorporated into our daily routines that we are not even aware of at any given moment. We truly are "given moments" and we should treasure them and make them work for us. Time to start anew.



*Camille Mayorga*

Sometimes we just need to update our "current" way of thinking, take a hard look at our "mind-chatter" and eliminate any negative thoughts we send ourselves. After all, our minds are listening when we tell ourselves that we can't remember names or where we put our car keys or that important piece of paper. Our mind is also listening when we tell ourselves we can never find a good parking place, or that we are always running late. Instead, fill your mind with positive thoughts and energy. Positive thoughts create around yourself an atmosphere favorable to the development of positive outcomes. "A man is what he thinks about all day long", famous words by Ralph Waldo Emerson, clearly states the power of the mind.

In keeping with fresh "current" thinking, the editorial staff here at *Currents* newsletter is committed to bringing you the most "current" in dental equipment news, software, products and events. We hope you enjoy this issue and are inspired to do some mental spring-cleaning.



By now you probably know how important the "first impression" really is. Your web site is your first impression on the Internet, make sure it's a good one.

YOU CAN BE SURE YOUR FIRST IMPRESSION IS A GREAT ONE WHEN CDS CREATES YOUR WEB SITE.

CDS is a full service new media company specializing in:

- ❖ Web Site Development - New or Redesign
- ❖ Graphic Design / Animation
- ❖ Corporate Identity
- ❖ Marketing - Print ads, Appearal, and much more
- ❖ Onsite Murals & Original Artwork
- ❖ Web Site Maintenance

"Their service is professional, fast and dependable. They encourage and coach us along in developing the content of our web site and marketing efforts. Our new web site has added to our bottom line"

Diagnostic Digital Imaging  
www.ddicentral.com

**We developed the new AADMRT.com website!**

**LET US SHOW YOU "THE POWER OF DESIGN"**

**PHONE: (916) 429 - 0694**



"feel the power of design"

e-mail: [Webmaster@CaDesignStudio.com](mailto:Webmaster@CaDesignStudio.com)

# News and Trends

## Other Dental Radiology Meetings

2004 American Association of Oral Maxillofacial Radiology Annual Session Nov. 3-7 Denver Colorado [www.aaomr.org](http://www.aaomr.org)  
 2004 European Congress of DentoMaxillofacial Radiology March 5-8, Malmo, Sweden [www.od.mah.se/eadmfr](http://www.od.mah.se/eadmfr)  
 2004 World Confrence of Oral Impantology, March 5-8, Honolulu, Hawaii, [www.wcoi6.com](http://www.wcoi6.com)  
 2005 International Congress of Dentomaxillofacial Radiology May 28th-June 1 Cape Town, South Africa, [www.iadmfr.org](http://www.iadmfr.org)



## New Feature on AADMRT.com

The public half of the ad rotator is complete and working on the AADMRT website ([www.aadmrt.com](http://www.aadmrt.com)). Instead of the long list of ads that has been accumulating over the past few months, the user is now presented with a random selection of three ads from the “pool” of advertising on the server. This selection changes with every page load (you can hit “refresh” a couple of times to see this). We can also prepare a report to the advertisers showing how many hits they have received through the AADMRT web site. The reports show day, month and yearly activity. This in turn can be very helpful for future and current advertisers with the AADMRT.

Ad		Views	Clicks	Click %
Aperio	Past Day:	21	0	0.00 %
	Past Week:	276	1	0.26 %
	Past Month:	652	7	1.07 %
	Past Year:	652	7	1.07 %
	Total:	652	7	1.07 %
CDS	Past Day:	31	0	0.00 %
	Past Week:	295	0	0.00 %
	Past Month:	760	2	0.26 %
	Past Year:	760	2	0.26 %
	Total:	760	2	0.26 %
Dolphin Imaging	Past Day:	31	0	0.00 %
	Past Week:	277	0	0.00 %
	Past Month:	730	3	0.41 %
	Past Year:	730	3	0.41 %
	Total:	730	3	0.41 %
Hitachi	Past Day:	26	0	0.00 %
	Past Week:	339	4	1.18 %
	Past Month:	739	5	0.69 %
	Past Year:	739	5	0.69 %
	Total:	739	5	0.69 %
Instrumentarium	Past Day:	34	0	0.00 %
	Past Week:	393	1	0.26 %
	Past Month:	762	2	0.26 %
	Past Year:	762	2	0.26 %
	Total:	762	2	0.26 %
Imaging Science	Past Day:	35	0	0.00 %
	Past Week:	391	1	0.26 %
	Past Month:	763	3	0.39 %
	Past Year:	763	3	0.39 %
	Total:	763	3	0.39 %
MME Consulting	Past Day:	35	0	0.00 %
	Past Week:	356	1	0.28 %
	Past Month:	742	9	1.21 %
	Past Year:	742	9	1.21 %
	Total:	742	9	1.21 %

## Spring Conference Nearly Sold Out

The 2004 spring meeting to be held in Costa Mesa, CA has very little seating left and we have outgrown our original room. We have worked with the hotel to have available to us a larger room so we can accomodate the high demands for this meeting, but if you plan on attending, you should sign up soon. This will be a sold out conference, and we will not be able to expand rooms again. Please see the insert in this publication for registration information, or go to [www.aadmrt.com](http://www.aadmrt.com) and click on future events.

## ORAD List

For those that want to discuss oral radiology, you may want to register on the Oral Radiology Discussion Group. This discussion group was started by Dr. Stuart White (Dental Radiologist, UCLA), and the topics include areas of interest in oral and maxillofacial radiology. If you are interested in subscribing, please log on to: <http://lists.ucla.edu/cgi-bin/mailman/listinfo/oradlist>



# News and Trends

## Imaging Prior to Implant Placement Suggested

Recent evidence has suggested that the use of radiographic images prior to implant placement can be beneficial, but until recently, there has been very little evidence supporting this. An article in the *Journal of the Massachusetts Dental Society* cites that a study showing that the use of tomography for pre-implant placement revealed an 87% agreement between the implant dimensions verified by the tomography, and the implant dimensions used in actual surgery. Without the use of tomography, there was only a 30% agreement between the implant dimensions, and those used during surgery. The results indicate that advanced imaging like tomography makes a positive difference in treatment planning for implants, it is cost effective, and has a significant dose reduction compared to spiral CT.

## Fall Conference- Salt Lake City, Utah

Make plans now to attend the annual fall conference this September 23rd-25th. The conference will be held in beautiful Salt Lake City Utah, surrounded by the picturesque Wasatch Mountains. This meeting will be held at the Park City Marriott located on historic Main Street. This conference will provide you with 12 CE credits and promises to be an interesting and informative session. Please see the advertisement for this meeting on page 8, or go to our website and click on future events for more specific details: [www.aadmrt.com](http://www.aadmrt.com)

## KODAK Digital Camera, Film Packets Earn Top *DentalTown* Awards

**ROCHESTER, N.Y.**, December 10—Two flagship products offered by Kodak's dental systems group have been voted by members of the *DentalTown* online community to win the 2003 Townie Choice Awards.

**KODAK INSIGHT** Intraoral Dental Film in SURESOFT Packets won the top award in the x-ray film category. In the extraoral camera category, the KODAK DX4900 Dental Digital Camera won the top award. Townie Choice Awards are selected by vote, open to the *DentalTown* online communities' 29,000 members.

### SURESOFT Packets

KODAK SURESOFT Packets feature a comfort edge that helps protect patients from the discomfort that can sometimes occur from the edge of a dental x-ray film packet. They are available with new KODAK INSIGHT Intraoral Dental Film, an F-speed film, which allows dental professionals to reduce radiation dosage for patients up to 60 percent compared with D-speed films.

### DX4900 Camera Kit

The KODAK DX4900 Dental Digital Camera Kit is designed to provide dentists with an easy, cost-effective way to create high-quality prints of dental images, including clinical views of teeth as well as patient portraits.

For more information on these or any other KODAK dental products, call 800-933-8031, visit Kodak's Health Imaging Dental Systems Group website at: [www.kodak.com/go/dental](http://www.kodak.com/go/dental).



# Spotlight Lab

Submitted by: *Betty Walls, R.T.*  
Wilson Radiographic Center of Austin  
Austin, Texas



I have been a Radiographic Technologist since 1967. I worked in hospital settings and family practice clinics for 10 years. Then in 1977 my love of the outdoors and whitewater canoe racing took me on a different path. I became the first River Outfitter in Texas. I led whitewater rafting trips in Texas, The Big Bend of Texas, and Mexico as well as renting canoes and rafts on the Guadalupe River. In 1987 I became involved in the politics of saving Texas Rivers and was instrumental in creating the first River District in Texas. In 1992 I sold my river companies and moved to Austin.



*Betty & Jeanette*

I assumed I would get a job with Texas Parks and Wild Life, but there were 800 applicants for one position so after 6 months of being told I was over qualified for most outdoor recreation jobs or environmental positions I looked at going back into the medical field. I was hired by a hospital that turned out to be the trauma center for 14 counties. So after being out of the field for 15 years I quickly got back up to speed. It didn't take long to realize trauma centers were not as exciting as they were when I was young. Luckily I saw an ad for Wilson Radiographic Center and wondered what do they do?

Wilson Radiographic Center of Austin was the first commercial dental lab to open in Texas.

Now there are four Wilson Labs in Houston, one in Dallas, and one in Plano, Texas. I was hired in 1992 and have gained much of my new skills in dental radiology by attending AADMRT conferences. We have a referral base of 40 orthodontist and some 60 general dentists. I act as chief technician and office manager and have cross-trained our receptionist Jeanette Hembrough to take radiographs and photos. She has completed the state dental x ray certification program and really enjoys the two roles she fills. Our team is small but we work great together. We have a dental assistant temp. who comes in to help during the summer when are schedule is the heaviest. We have the Instrumentarium Pan Machine and a Quint Sectograph, old but still the workhorse of the industry in my opinion. We also utilize the Dolphin Imaging system for our intra-oral photos.

I really enjoy this job and the interaction with children in particular. One of my most memorable moments was when I was working on a child and asked her to smile for the photos. She was only 8 years old but she said "my mom told me not to smile with my teeth showing because they are ugly." I felt so sorry that she felt that way but I knew when she finished her orthodontic treatment she would have a great smile. Knowing that we contribute to having that happen makes me feel really good. I also appreciate the support and training that AADMRT provides and how the other labs share their knowledge with us when we attend conferences. And for fun I still go whitewater rafting and saltwater fishing.

**Table 1**

**Comparison of Maxillofacial CBYT Devices**

Parameters	Newtom 9000	Newtom Plus	3DX Accutomo	ISI/CAT	Hitachi MercuRay	MedicalCT
X-ray beam	cone	cone	cone	cone	cone	fan
Sensor detector	area image intensifier CCD	area image intensifier CCD	area image intensifier CCD	area amorphous silicon flat-panel detector	area image intensifier CCD	linear Solid state or gas
Grayscale	8 bit	12 bit	8 bit	12 bit	12 bit	12 bit
Voxel size (mm <sup>3</sup> )	0.265	0.07-0.20 [variable]	0.125	0.4 typical 0.2 minimum	0.1	0.3 <sup>16</sup>
X-ray source			direct current	direct current		
anode	fixed	fixed	fixed	fixed	fixed	rotating
kVp range	110 [fixed]	110 [max]	60-80	120	70-100	110-140
mA	10-15 max		1-10	1-3	to 15 max	80-300
Image acquisition	panoramic type single 360° rotations	panoramic type single 360° rotations	panoramic type single 360° rotations	panoramic type single 360° rotations	panoramic type single 360° rotations	axial slices multiple 360° rotations
Patient position	supine	supine	seated	seated	seated	supine
Image area	maxillofacial	maxillofacial	maxillofacial	maxillofacial	maxillofacial	entire body
Dimensions (cm)	13x13 (height x diameter)	22x25 (height x diameter)	3.0x4.0 (height x diameter)	2 versions 11x17 17x17 (height x diameter)	D mode 5.12 P mode 11.7 C mode 15.0 (height only)	varies/exam
Imaging session (sec)	75	30 or less	17	40 or less	9.6	varies/exam
Effective dose* (max & mand) mSv	.04-.05 <sup>15</sup>	.01-.02	.0074	not available	not available	.289 low <sup>17</sup> .723 high <sup>17</sup>
Auto exposure control	smart scan yes	smart scan yes	no	not reported	not reported	no
Commercially available	yes	end of 2003/early 2004	yes	projected 2004	not reported	yes

\*Effective dose as reported is provided from the manufacturer and is not truly comparable to other devices due to variable such as different size of imaging volumes and operational settings of the devices.

### Principles and Brief History of X-ray Computed Tomography

CT scanners consist of an X-ray source and detector mounted on a rotating gantry (Fig. 2). During one rotation of the gantry, the detector detects the flux,  $I$ , of X-rays that have passed through the patient. The attenuation of monochromatic X-rays in homogenous objects is governed by:  $I = I_0 \exp(-\mu x)$ , where  $I_0$  is the X-ray intensity without the object,  $x$  is the length of the X-ray path through the object, and  $\mu$  is the linear attenuation coefficient of the material at the X-ray energy employed. For inhomogeneous objects, like the human body, the attenuation of X-rays consequently can be described  $I = I_0 \exp(-\int \mu(x) dx)$ . By taking the logarithm of the fluence  $-\log(I/I_0)$  one obtains line integrals of the linear attenuation coefficients. These integrals constitute so-called “raw data” that are then fed into an image reconstruction method that generates cross-sectional images whose pixel values correspond to linear attenuation coefficients. The theoretical background for tomographic image reconstruction was laid out in 1917 when Radon established that a 3D object can be reconstructed from an infinite set of 2D projections obtained at varying angles around the object.

The resulting attenuation coefficients are usually expressed with reference to water, and are given in Hounsfield units (HU):

$$HU_{\text{patient}} = 1000 \times \frac{\mu_{\text{patient}} - \mu_{\text{water}}}{\mu_{\text{water}}}$$

The first CT scanner was developed in 1967 by Sir. Godfrey N. Hounsfield, an engineer at EMI. Since then, CT technology rapidly underwent four developmental generations. The first generation of CT scanners used a single detector element to capture a beam of X-rays, corresponding to the integral of linear attenuation coefficients along a single line. It then translated horizontally to acquire the next line integral. After acquiring all the line integrals for a given position of the X-ray source, both the detector and source rotated one degree—a design known as the “translate-rotate” or “pencil-beam” scanner. Hounsfield’s unit belonged to this generation, as did the first commercial CT scanners introduced in 1972. Interestingly, these first generation CT scanners were designed to scan the head only.

A second generation of CT systems was introduced in 1975. These systems, also known as “hybrid” machines, used more than one detector and used small fan-beam, as opposed to pencil-beam, scanning. Like the first generation of CT scanners, these scanners also used a translate-rotate design, and were for the most part head-only scanners. While the first iterations of full body CT scanners also incorporated the translate-rotate design, image quality was poor due to patient motion artifacts caused by the significant amount of time required to take the scan.



Figure 1: MiniCAT/CAT scanner developed by Xoran Technologies, Inc. and Imaging Sciences, Inc.




Figure 2: Principles of conventional (“fan beam”) CT (left) and cone beam CT (right)

You've updated  
your equipment.  
Isn't it time to update  
your delivery?




Securely Share: Digital Photos Radiographic Images Cephalometric Tracings\*

**Introducing Dolphin Imaging Version 9, featuring AnywhereDolphin™**—the Internet-integrated imaging software that is revolutionizing the way you practice and deliver records.

Using the same data encryption technology as banking institution sites, **Version 9** lets you store, organize and share records through an encrypted, password-protected online account.

For more information on how **Version 9** can enhance the way you practice, visit [www.dolphinimaging.com](http://www.dolphinimaging.com) or call us at 1-800-548-7241.

\* Available Spring 2004 © 2003 Dolphin Imaging & Management Solutions



Third generation CT scanners appeared in 1976 and are the systems most widely used today. These scanners use a large, arc-shaped detector that acquires an entire projection without the need for translation. This rotate-only design, frequently referred to as “fan-beam”, utilizes the power of the X-ray tube much more efficiently than the previous generations.

Fourth generation scanners shortly followed third generation scanners, replacing the arc-shaped detector with an entire circle of detectors. In this design the X-ray tube rotates around the patient, while the detector stays stationary. Since these fourth generation scanners tend to be more expensive and suffer from higher levels of Compton scatter artifacts, most of the commercially available CT scanners today are third generation scanners.

After an initial period of rapid development, CT technology quickly became mature, and it was not until the early 1990s that CT research began anew. Recent advances in CT include multirow detectors and spiral scanning. Multirow scanning allows for the acquisition of several cross-sectional slices at the same time, reducing scanning times. Today’s state-of-the-art scanners have 64 rows of detectors. Spiral (helical) scanning incorporates a moving table with the rotating X-ray tube, with the net effect that the X-ray tube describes a helical path around the patient.

Conventional CT scanners are large and expensive systems designed for full-body scanning at a high speed to minimize artifacts caused by movement of the heart, lungs, and bowels. They are not well-suited for dedicated dentomaxillofacial imaging, where cost considerations are important, space is often at a premium, and scanning requirements are limited to the head. The advent of cone beam CT (CBCT) technology has paved the way for the development of relatively small and inexpensive CT scanners dedicated for use in dentomaxillofacial imaging.

### **Principles and Brief History of Cone Beam Computed Tomography**

Cone-beam CT scanners utilize a two dimensional detector (Fig. 2), which allows for a single rotation of the gantry to generate a scan of the entire region of interest, as compared to conventional CT scanners whose multiple “slices” must be stacked to obtain a complete image. In comparison with conventional fan-beam or spiral-scan geometries, cone-beam geometry has higher efficiency in X-ray use, inherent quickness in volumetric data acquisition, and potential for reducing the cost of CT. Conventional fan-beam scans are obtained by illuminating an object with a narrow, fan-shaped, beam of X-rays. The X-ray beam generated by the tube is focused to a fan-shaped beam by rejecting the photons outside the fan, resulting in a highly inefficient use of the X-ray photons. Further, the fan-beam approach requires reconstructing the object slice-by-slice and then stacking the slices to obtain a 3D representation of the object. Each individual slice requires a separate scan and separate 2D reconstruction. The cone beam technique, on the other hand, requires only a single scan to capture the entire object with a cone of X-rays. Thus, the time required to acquire a single cone-beam projection is the same as that required by a single fan-beam projection. But since it takes several fan beam scans to complete the imaging of a single object, the acquisition time for the fan beam tends to be much longer than with the cone beam. Although it may be possible to reduce the acquisition time of the fan beam method by using a higher power X-ray tube, this increases the cost and size of the scanner as well as the electric power consumption, thus making the design unsuitable for a compact scanner.

A number of groups have worked on developing task-specific CBCT scanners over the past two decades. Computed tomography angiography (CTA), in particular, has been an active area of investigation due to its lenient requirements for contrast resolution and strict requirements on spatial resolution - a natural fit for CBCT. The first CBCT scanner ever to be built was built for angiography among other tasks at Mayo in 1982 [1].

# MME will help you apply technology to make your lab more profitable!

There are many amazing tools available today that can **add value and profit** to your dental x-ray lab. MME is **experienced** with these tools and applying them to make an **efficient and successful** dental x-ray lab.

## MME can help you:

- ✓ Optimize communication by connecting your office to email and the Internet
- ✓ Share office resources by networking lab computers
- ✓ Integrate Digital Photography solutions through digital cameras, photographic printers, imaging software, data storage and backup
- ✓ Incorporate Digital X-ray solutions
- ✓ Invest in an Internet & Community presence by adding a website to your business

Call us today or find us online!

Phone 916-419-1102  
www.mmeconsulting.com

We put the  
"digital"  
in ddi labs.



12 McKilt Ct, Sacramento, CA 95835-1334 email: information@mmeconsulting.com

## AADMRT COMMITTEE REPORTS

The AADMRT Board of Directors held it's meeting in Sacramento CA, on February 7<sup>th</sup> 2003; Merry Hampton, Executive Secretary, submits the general membership minutes.

The Board welcomes two new members: Kevin Fox, who will take over as Recruitment Chair, and Tracey Saucier, will take over the By-Laws. Craig Dial has agreed to stay on as President. Duane Perry will continue as Vice President.



From left to right: Jeannie Herriott, Tracey Saucier, Kathleen Cox, Duane Perry, Bart Webb, Kevin Fox, Craig Dial

**Membership:** 184 current members, 49 non-paid .New members in 2003: 53,

New members 2004: 4 , Labs in database:153

**Nominating:** needs three new members for 2005.

**Website:** Paypal is now available for credit card payments online for membership dues.

**Bi-Laws:** Board agreed to continue using the term "Technicians" vs. changing to "Technologist".

**Continuing Education:** currently have 7 students enrolled in our school.

**2004 Fall Conference:** Park City, Utah, September 23-25. Duane Perry Chairperson.

**2005 Conference Sites:** Spring -San Fransisco, CA; Fall -Phoenix, AZ

**New Business:** board agreed to place an advertisement in one or more of the National Dental Journals.

Fahrig et al [2] [3] have developed a CBCT system based on an image intensifier and C-arm for use in angiography. Wiesent [4] have also developed a C-arm plus image intensifier system for interventional angiography. Saint-Felix et al [5] developed a CTA CBCT system based on the gantry of a conventional CT scanner which reconstructs vasculature from a set of digitally subtracted angiography (DSA) images. Ning et al. have developed a CBCT angiography imager based on GE 8800 CT scanner with an image intensifier - CCD chain and later with a flat-panel detector [6] [7] [8]. Schueler et al. have developed a CBCT CTA scanner based on a biplanar C-arm system [9]. Kawata et al. also developed a CBCT CTA system [10].

Jaffray and Siewerdsen have developed a CBCT system for radiotherapy guidance based on an amorphous silicon (a-Si:H) flat-panel detector [11] [12, 13]. Cho et al. have also developed a CBCT system for radiotherapy applications [14]. Efforts are being made towards dedicated CBCT-based imaging systems for mammography [15].

Although CBCT has existed for over two decades, its true potential has not yet been fully tapped. Only recently has it become possible to develop CBCT clinical systems that are both inexpensive and small enough to be used in OR, medical offices, emergency rooms, and intensive care. Four technological and application-specific factors have converged to make this possible. First, compact and high-quality flat-panel detector arrays were developed. Second, the computer power necessary for cone-beam image reconstruction has become widely available and is relatively inexpensive. Third, x-ray tubes necessary for cone-beam scanning are orders-of-magnitude less expensive than those required for conventional CT. Fourth, by focusing on head/neck scanning only, one can eliminate the need for sub-second gantry rotation speeds that are needed for cardiac and thoracic imaging. This significantly reduces the complexity and cost of the gantry.

In short, cone beam CT is ideally suited for high quality and affordable CT scanning of the head and neck in dentomaxillofacial applications.




Figure 3: Axial, coronal, and sagittal slices of a maxillofacial scan obtained with MiniCAT™ / I-CAT.

## Pacific Dental Mounts & Supplies



- ◆ X-Ray and Ortho Photo Mounts
- ◆ Polyester Sleeves and Sheets, sturdy .5 mil thick, cut to almost any size
- ◆ Mailing and Case Storage Envelopes, thumbcut to almost any size
- ◆ Tracing Acetate
- ◆ Before and After Folders, standard or custom foil-stamped
- ◆ Labels

CALL FOR CATALOG OR SAMPLES

800-323-7068 707-254-0698

Napa, California

This value of using cone beam CT for dedicated dentomaxillofacial imaging has been recognized by a number of researchers, and several commercial systems are available commercially from Quantitative Radiology ("NewTom 9000" [3]), Hitachi, J. Morita Co. [4,5], as well as from the collaboration of Xoran Technologies and Imaging sciences International (MiniCAT™ / I-CAT (formerly DentoCAT™) [6])

Continued on page 16

# CB MercurRay™

Changing the face of dental imaging

Hitachi has taken its expertise in conebeam and CT technologies and applied it to the field of dentistry with the CB MercurRay. The result is an x-ray scanner with a 360° rotational sequence, enabling 2D images and 3D volume reconstructions for excellent planning and diagnostic support.

## Advanced Digital Imaging Specialized for Dentistry

- Distortion free panoramic x-rays
- TMJ linear digital tomograph series
- Cephalometric x-rays including soft tissue
- Isotropic CT volume (custom slicing in any direction)
- 3D interactive volume rendering of the patient anatomy (skulls and jaw)
- Paranasal sinuses and airway visualization

## Features for Superior Imaging and Ease of Use

- Simple acquisition with user-friendly software
- High resolution
- Reduction of artifacts
- Distortion free
- Excellent contrast resolution (for bone and soft tissue)
- True quantitative data (with isotropic resolution)



## Fast and Comfortable for Patient Convenience

- Less than 10 second scan time
- Patient in normal seated position
- Quick system set-up
- Adjusts to patient size
- Preview mode to ensure correct positioning

All provided with Hitachi's renowned  
service, applications and marketing support

# HITACHI

HITACHI MEDICAL SYSTEMS AMERICA, INC.

1959 Summit Commerce Park • Twinsburg, Ohio 44087 • 1-800-830-3106 • [www.hitachimed.com](http://www.hitachimed.com)

### Flat-Panel Detector based on Amorphous Silicon

The quality of the reconstructed CT image depends significantly on the quality of the data that is acquired during the scan. The characteristics of the two-dimensional detector used by a cone beam CT scanner affect therefore directly the quality of the CT scan.

MiniCAT™ scanner deploys a flat-panel detector constructed from hydrogenated amorphous silicon (aSi:H). These detectors are essentially self-scanned arrays of n-i-p photodiodes and thin-film transistors switches with a scintillator layer. Essentially the same underlying technology is used to construct flat-panel computer monitors and large-area document imagers. Flat-panel detectors are replacing film and image-intensifiers in conventional radiography and fluoroscopy applications.

The characteristics of aSi:H based flat panel detector arrays make them a much better choice for cone beam imaging than an alternative solution: image intensifier coupled with a CCD. Image intensifiers create geometric distortions that must be addressed when processing the data later in the software, while flat-panel detectors do not suffer from this problem. Also, flat panel detectors afford a greater dynamic range than that offered by the image intensifier + CCD camera approach.

### Sample Images

Figures 3-6 illustrate the high spatial resolution that can be obtained with CBCT (MiniCAT™/I-CAT) and compare it to that of conventional CT scanner.




Figure 4 Comparison of a conventional CT scan (left) with a CBCT scan (right). Note the higher resolution along the z-axis for the CBCT scan.



Figure 6 3D renderings of a MiniCAT scan of an image-guided maxillary surgery phantom. The high isotropic spatial resolution makes it ideal for example for image-guided functional endoscopic sinus surgery.

### Future of Cone Beam CT

Its high spatial resolution, smaller size and lower cost has made CBCT a natural fit dentomaxillofacial imaging. Those same attributes will drive the adoption of CT in other markets as well. Xoran Technologies is successfully marketing the MiniCAT™ scanner to otolaryngology offices for in-house imaging of the sinuses and temporal bone. In the future, integration of CBCT with image guided surgery systems for functional endoscopic sinus surgery is anticipated (a scan of an image guided surgery phantom is shown in Figure 6.) 3D intraoperative imaging devices based on CBCT are also in development.

At the other end of the spectrum, the ‘big four’ companies that manufacture conventional CT scanners are working on mating cone beam technique with the scanners from their product lines. Once accomplished, this would allow the physicians to ‘freeze’ the motion of the heart for cardiac imaging. The latest generation of conventional CT scanners has 64 slices and can cover 4 cm axially. It is anticipated that the need for sophisticated cardiac imaging tools will drive the introduction of truly cone beam helical scanners within several years.

conclusion on page 17

As CBCT scanning is finding more and more commercial applications in medicine, dentomaxillofacial radiology stands as the privileged field that has driven growth of this exciting technology out of the R&D infancy into the commercial maturity.




Figure 5 3D renderings of a high-resolution isotropic (voxel size 0.2 mm) scan of a dry mandible.

**Bibliography**

1. Robb, R.A., *Dynamic Spatial Reconstructor: An X-ray Video Fluoroscopic CT scanner for dynamic volume imaging of moving organs.* IEEE Trans Med Im, 1982. **MI-1**(1): p. 22-23.
2. Fahrig, R., et al., *Use of a C-arm system to generate true three-dimensional computed rotational angiograms: preliminary in vitro and in vivo results.* American Journal of Neuroradiology, 1997. **18**: p. 1507-14.
3. Fahrig, R., et al. *Computed rotational angiography: system performance assessment using in vitro and in vivo models.* in *SPIE Conference on Medical Imaging.* 1998. San Diego: SPIE.
4. Wiesent, K., et al., *Enhanced 3D reconstruction algorithm for c-arm systems suitable for interventional procedures.* IEEE Trans Med Im, 2000. **19**(5): p. 391-403.
5. Saint-Felix, D., et al., *In vivo evaluation of a new system for 3D computerized angiography.* Phys. Med. Biol., 1994. **39**: p. 584-95.
6. Ning, R., et al., *Flat panel detector-based cone beam volume CT angiography imaging: system evaluation.* IEEE Trans Med Im, 2000. **19**(9): p. 949-63.
7. Wang, X., *Volume Tomographic Angiography.* 1997, University of Rochester.
8. Ning, R., et al. *Real time flat panel detector-based volume tomographic angiography imaging: detector evaluation.* in *SPIE Conference on Medical Imaging.* 2000. San Diego: SPIE.
9. Schueler, B.A., et al., *Three-dimensional Vascular Reconstruction with a Clinical X-ray Angiography System.* Acad Radiol, 1997. **4**(10): p. 693-699.
10. Kawata, Y., N. Niki, and T. Kumazaki, *Measurement of blood vessel characteristics for disease detection based on cone beam CT images.* IEEE Trans Nuc Sci, 1996. **43**(6): p. 3348-54.
11. Jaffray, A.D. and J.H. Siewerdsen, *Cone-beam computed tomography with a flat-panel imager: Initial performance characterization.* Medical Physics, 2000. **27**: p. 1311-1323.
12. Siewerdsen, J.H. and D.A. Jaffray, *Cone-beam computed tomography with a flat-panel imager: magnitude and effects of x-ray scatter.* Medical Physics, 2001. **28**: p. 220-31.
13. Siewerdsen, J.H. and D.A. Jaffray, *Cone-beam computed tomography with a flat-panel imager: effects.* Medical Physics, 1999. **26**: p. 2635-47.
14. Cho, P.S., R.H. Johnson, and T.W. Griffin, *Cone-beam CT for radiotherapy applications.* Phys. Med. Biol., 1995. **40**: p. 1863-83.
15. Ning, R. and B. Chen. *Cone beam volume CT mammographic imaging: feasibility study.* in *Presented at Proceedings of SPIE - The International Society for Optical Engineering.* 2001.

**We're Moving!** Same Great Service!  
Same Great Quality!

**P&G PRODUCTS**  
*P&G*  
Custom Design

**P&G Products**  
Toll Free (800) 367-8847  
(805) 226-7333 • (805) 226-7337 Fax


**AADMRT Supporter for 28 years!**

- X-Ray Mounts • Photo Mounts • Polyester Sleeves
- Polypropylene Sleeves
- Tracing Acetate • Case Envelopes
- Digital, Presentation and Before & After Folders

*New Address - P.O. Box 3680 / Paso Robles, CA 93447*

**FULLY INTEGRATED PATIENT SCHEDULING**

**Quick & Easy!**



- ▶ One-time patient data entry & schedule
- ▶ Convenient drag & drop
- ▶ Easy follow-up/rescheduling
- ▶ Print daily schedule & notes
- ▶ Optimize equipment usage

**Radiology Lab Management System**

- Network ready
- NEW! Daily notepad
- Billing and accounts receivable for patients, doctors, insurance companies - always updated
- Daily and periodic summaries
- X-ray cards, photo labels, receipts
- Superbills
- Schedules and tracks deliveries
- Optional automated e-mail delivery for digital photos

*"RLMS is fast, efficient and easy to learn."  
Neville Jones, Deputy Manager (Lab Eng)*

*"Finally, someone realized the needs and developed a quality solution to our problem. The savings started immediately."  
Ken J. Technick, Director (Lab Eng)*



**DESIGN TECHNOLOGY**  
Information and Technology Resources at Work for You.  
**(303) 756-4333**  
www.designtech.net

# Lab Products

## Dental Radiography Laboratory Manual

By Sandra Slack Olson, RDH, MEd, CAGS, This manual provides practical material specific to the dental radiography laboratory experience. Each chapter contains learning objectives, review exercises, and learning activities, enabling the reader to further their competency level. 251 pp. ISBN: 0721664555 Price: \$34.95 [www.quintpub.com](http://www.quintpub.com) 1-800-621-0387

## Wireless CDR

Wireless active pixel sensors designed to give you everything you expect from a sensor like high quality images, reduced radiation, and simple to use plus no cord! For more information go to: [www.schichteck.com](http://www.schichteck.com)




## Fixer Treatment Systems

Compact chemical processing unit keeps office in compliance with federal, state and local EPA standards. Easy to use, can be set up wherever there is a sink or drain. [www.smartpractice.com](http://www.smartpractice.com) Price: \$166.60



## Replacement Parts and Service

For Dent-X/Philips, Air Techniques and Hope/Litton/Allied Automatic X-Ray Film Processors. Replacement parts include design, material and manufacturing enhancements to improve the quality and reliability of your X-ray film processing. [www.smithdental.com](http://www.smithdental.com)



## HIPAA Forms

The dental record can help with your HIPAA compliance. We have the forms you need, and free copies for your review. Just call:

1-800-243-4675, or on line at: [www.dentalrecord.com](http://www.dentalrecord.com)



## Sterilization Pouches

Crosstex sterilization pouches are ADA accepted. Our pouch features an extra-thick, fluid-resistant paper stock utilizing a triple seal, PLUS a fourth "Secure Seal" at the chevron for thorough protection of contents. Pouches have easy-to-read color indicators for steam, chemical vapor, or gas sterilization, which assures sterile conditions have been achieved. Pouch Self Seal 12x18 \$39.99 [www.henryschein.com](http://www.henryschein.com)

## Sell or Buy used X-Ray Equipment





DentalQuip, International specializes in the export of new and refurbished dental equipment in Mexico and other Latin America Countries.

DentalQuip, International is the sole representative for Forest Medical products in Mexico. Forest Medical Products is a leading dental equipment manufacturer in the U.S.A.

[www.dentalquip.com](http://www.dentalquip.com)





### Newsletter Deadlines

Winter Issue Deadline: December 1  
Spring Issue Deadline: March 1  
Summer Issue Deadline: June 1  
Fall Issue Deadline: September 1

For Advertising Information :  
Advertising@aadmrt.com  
For Article Submission info:  
Editor@aadmrt.com

Or Visit our Website:  
[www.aadmrt.com](http://www.aadmrt.com)

### Future AADMRT Event Calender

**2004 Spring Meeting**  
Costa Mesa, CA  
March 27

**2004 Fall Conference**  
Salt Lake City, UT  
September 23-25

**2005 Spring Meeting**  
San Francisco, CA  
date- TBA

**2005 Fall Conference**  
Phoenix, AZ  
date- TBA

*Check web site for more specific details*

*The Organization acknowledges  
the sponsors of this Newsletter  
with thanks:*

*Hitachi Medical Systems  
Aperio Services  
Pacific Dental Mounts  
MME Consulting  
California Design Studio  
Panoramic Dental  
Imaging Sciences International  
Design Technology  
P&G Products  
Dolphin Imaging  
Kodak*

**[www.aadmrt.com](http://www.aadmrt.com)**

### X-Ray Lab & Imaging Currents

AADMRT  
1 Scripps Drive #101  
Sacramento, CA 95825

Place Stamp  
Here

